

## **DENTAL UNIT** **MODEL: MSQDNT1**



### **Standard Equipment:**

- Rotary Tempered glass spittoon x1 set
- Auto spitton flush and cup filler control system x1 set
- Water tank (air pressurized) (built –in) x1 set
- Operating light x1 set
- Instrument tray with air break x1 set
- Saliva ejector with high and low suction x1 set
- 3 way syringe (cold/hot water) x2 set
- 24V noiseless DC motor chair x1 set
- Doctor stool x1 set
- Rotary ceramic spittoon x1 set
- 

### **CENTURION**

Stand 566B  
6661 Rooikeefflap Close,  
Centurion, Pretoria  
P.O. Box 19, Brits, 0268  
Telephone (083) 263 7213  
Fax 086 574 0637  
E-mail: info@msquaredmedical.co.za

**Options:**

- **Imported 4 hole high speed air turbine hand piece**
- **Imported 4 hole low speed air motor hand piece**
- **Scaler**
- **LED curing light**
- **Endoscope system**
- **Leather padded cushion seat**

**Specifications:**

- power voltage: 220V / 50HZ
- motor voltage: 24V
- water pressure: 0.2MPa-0.4MPa
- air pressure: 0.5MPa-0.8MPa



**CENTURION**

Stand 566B  
6661 Rooikeeflap Close,  
Centurion, Pretoria  
P.O. Box 19, Brits, 0268  
Telephone (083) 263 7213  
Fax 086 574 0637  
E-mail: info@msquaredmedical.co.za