

ELECTROSURGICAL UNIT SURTRON 50 D



DENTAL
DERMATOLOGY
AMBULATORY SURGERY
FIRST AID
VETERINARY
VASCULAR SURGERY



SURTRON® 50 D is a radiofrequency electrosurgical suitable for small monopolar surgery.

SURTRON® 50 D through the functions selection makes possible pure cut, coagulation cut, and pure coagulation. The digital display of power erogation allows safety and very reproducible applications.

SURTRON® 50 D permits use a professional surgery thanks to the stability functions and the safety system adopted.

CENTURION

Stand 566B
6661 Rooikeyflap Close,
Centurion, Pretoria
P.O. Box 19, Brits, 0268
Telephone (083) 263 7213
Fax 086 574 0637
E-mail: info@msquaredmedical.co.za

Main Features:

SURTRON® 50 D	
Maximum output power CUT (W)	50 → 400 Ω
Maximum output power CUT/COAG 1 (W)	45 → 400 Ω
Maximum output power CUT/COAG 2 (W)	40 → 400 Ω
Maximum output power COAG (W)	40 → 400 Ω
Working frequency	600 kHz
Mains voltage (Vac)	115-230
Mains frequency Hz)	50-60
Electrical class (EN60601-1)	I CF
Neutral electrode	-F-
Weight Kg	2,8
Size HxLxD mm	180x100x180

Standard Accessories

SURTRON® 50 D	
Unit electricalsurgery code	10100.05A
Cable for rod neutral electrode	•/1
Kit of assorted electrode length 5cm (10pcs)	•/1
Manual of instructions	•/1
Not water-proof foot switch	•/1
Power supply cable 3x1mm SIE-IEC – 2mt	•/1

CENTURION

Stand 566B
6661 Rooikeyflap Close,
Centurion, Pretoria
P.O. Box 19, Brits, 0268
Telephone (083) 263 7213
Fax 086 574 0637
E-mail: info@msquaredmedical.co.za